


**OLYMPUS**

Your Vision, Our Future

Endo**Therapy**

# EZSHOT 3 PLUS


Access Difficult Locations with Ease and Confidence



## UNCOMPROMISED ACCESS TO THE PANCREATIC HEAD


### Needle Flexibility Provides Unparalleled Access

- Reduced insertion force even when scope is angulated
- No need for scope position adjustments to compensate for larger diameter needle
- Improved access to lesions in difficult locations





### Access from Wider Angles

Deep tissue can be punctured regardless of scope angle.




### Resilient Needle Shape

Nitinol retains its shape without any deformation even after multiple passes. This helps ensure the needle follows a predictable and consistent trajectory.



### Side Hole Option

Both models are also available with a side hole which helps to collect the sample within the needle.



# EZSHOT 3 PLUS

## Specifications

Article Name	Order Number	Quantity	Min. Working Channel Ø	Working Length	Max. Needle Length	Needle Ø	Stylet Tip Shape
NA-U200H-8019	N5409330	5	2.8 mm	1400 mm	80 mm	19G	Rounded
NA-U200H-8019S	N5409530	5	2.8 mm	1400 mm	80 mm	19G with side hole	Rounded
NA-U200H-8022	N5409730	5	2.8 mm	1400 mm	80 mm	22G	Rounded
NA-U200H-8022S	N5409930	5	2.8 mm	1400 mm	80 mm	22G with side hole	Rounded

## Standard Set Contents

A set: 5 pcs/box



Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

**OLYMPUS**

**OLYMPUS EUROPA SE & CO. KG**  
Postbox 10 49 08, 20034 Hamburg, Germany  
Wendenstrasse 14-18, 20097 Hamburg, Germany  
Phone: +49 40 23773-0, Fax: +49 40 233765  
www.olympus-europa.com